

PPP1R15B Rabbit pAb


ABZ1528

Product Overview


Name	PPP1R15B Rabbit pAb
Accession(Primary)	Q5SWA1
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 79 kDa; Observed MW: 79 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:10000

Target information (Q5SWA1)


Gene Name	PPP1R15B
About Gene Symbol	




Western blot analysis of PPP1R15B in NIH3T3 lysates using PPP1R15B antibody.



Immunohistochemistry analysis of paraffin-embedded Human liver using PPP1R15B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human stomach using PPP1R15B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human stomach using PPP1R15B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.