

CD62P Rabbit pAb


ABZ1529

Product Overview


Name	CD62P Rabbit pAb
Accession(Primary)	P16109
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 91 kDa; Observed MW: 91 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:10000

Target information (P16109)


Gene Name	SELP
About Gene Symbol	




Western blot analysis of CD62P in NIH3T3, KB lysates using CD62P antibody.



Western blot analysis of CD62P in 3T3 lysates using P-Selectin antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsils using CD62P antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human tonsils using CD62P antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.