

Growth Hormone Rabbit pAb

ABZ1532


Product Overview

Name	Growth Hormone Rabbit pAb
Accession(Primary)	P01241
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:50-1:200; ELISA-1:10000-1:20000; IF-1:50-1:200


Target information (P01241)

Gene Name GH1/2


About Gene Symbol




Immunohistochemistry analysis of paraffin-embedded Human placenta using Growth Hormone antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using Growth Hormone antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human brain using Growth Hormone antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human brain using Growth Hormone antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.