


IL-37 Rabbit pAb

ABZ1533


Product Overview

Name	IL-37 Rabbit pAb
Accession(Primary)	Q9NZH6
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human
Recommended Dilution	IHC-1:50-1:200; ELISA-1:10000-1:20000; IF-1:50-1:200


Target information (Q9NZH6)

Gene Name IL37**About Gene Symbol**


Immunohistochemistry analysis of paraffin-embedded Human small intestine using IL-37 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using IL-37 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human coloncancer using IL-37 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human coloncancer using IL-37 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.