

PREP Rabbit pAb


ABZ1535

Product Overview


Name	PREP Rabbit pAb
Accession(Primary)	Q8N6D4
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 81 kDa; Observed MW: 81 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (Q8N6D4)


Gene Name	PREP
About Gene Symbol	



Western blot analysis of PREP in Jurkat lysates using PREP antibody.



Immunohistochemistry analysis of paraffin-embedded rat brain using PREP antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded mouse brain using PREP antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.