

Transferrin Rabbit pAb

ABZ1540

Product Overview


Name	Transferrin Rabbit pAb
Accession(Primary)	P02787
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 77 kDa; Observed MW: 77 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (P02787)


Gene Name	TF
About Gene Symbol	



Western blot analysis of Transferrin in K562 lysates using Transferrin antibody.



Western blot analysis of Transferrin in HELA lysates using TF antibody.



Immunohistochemistry analysis of paraffin-embedded rat liver using Transferrin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.