

Neuropilin 1 Rabbit pAb


ABZ1543

Product Overview


Name	Neuropilin 1 Rabbit pAb
Accession(Primary)	O14786
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 103 kDa; Observed MW: 100 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (O14786)


Gene Name	NRP1
About Gene Symbol	



Western blot analysis of Neuropilin 1 in 293 lysates using Neuropilin 1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human kidney using Neuropilin 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human kidney using Neuropilin 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.