

Integrin alpha 4 Rabbit pAb


ABZ1544

Product Overview


Name	Integrin alpha 4 Rabbit pAb
Accession(Primary)	Q13349
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 127 kDa; Observed MW: 127 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (Q13349)


Gene Name	ITGAD
About Gene Symbol	



Western blot analysis of Integrin alpha 4 in 4T1 293 lysates using Integrin alpha 4 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung using Integrin alpha 4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human spleen using Integrin alpha 4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.