

HRP3 Rabbit pAb


ABZ1545

Product Overview


Name	HRP3 Rabbit pAb
Accession(Primary)	Q9Y3E1
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 23 kDa; Observed MW: 23 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (Q9Y3E1)


Gene Name	HDGFL3
About Gene Symbol	



Western blot analysis of HRP3 in NIH3T3, PC-12 lysates using HRP3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human Lymphnodes using HRP3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human Lymphnodes using HRP3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.