

RNA Polymerase III Subunit C7 Rabbit pAb

ABZ1546

Product Overview


Name	RNA Polymerase III Subunit C7 Rabbit pAb
Accession(Primary)	O15318
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 26 kDa; Observed MW: 26 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (O15318)


Gene Name	POLR3G
About Gene Symbol	



Western blot analysis of RNA Polymerase III Subunit C7 in PC-12 lysates using RNA Polymerase III Subunit C7 antibody.



Immunohistochemistry analysis of paraffin-embedded mouse brain using RNA Polymerase III Subunit C7 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded mouse brain using RNA Polymerase III Subunit C7 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.