

NOL5A Rabbit pAb


ABZ1552

Product Overview


Name	NOL5A Rabbit pAb
Accession(Primary)	O00567
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 66 kDa; Observed MW: 66 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:10000

Target information (O00567)


Gene Name	NOP56
About Gene Symbol	



Western blot analysis of NOP56 in HEPG2 , Hela lysates using NOP56 antibody.



Immunohistochemistry analysis of paraffin-embedded Human placenta using NOP56 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human placenta using NOP56 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.