

CD16 Rabbit pAb

ABZ1554


Product Overview

Name	CD16 Rabbit pAb
Accession(Primary)	P08637
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Rat
Recommended Dilution	IHC-1:50-1:200; ELISA-1:10000-1:20000; IF-1:50-1:200


Target information (P08637)

Gene Name FCGR3A/FCGR3B


About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded Human kidney using CD16 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using CD16 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human brain using CD16 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.