


IL18R Beta Rabbit pAb

ABZ1555


Product Overview

Name	IL18R Beta Rabbit pAb
Accession(Primary)	O95256
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Recommended Dilution	IHC-1:50-1:200; ELISA-1:10000-1:20000; IF-1:50-1:200


Target information (O95256)

Gene Name IL18RAP**About Gene Symbol**

Immunohistochemistry analysis of paraffin-embedded Human kidney using IL-18 Receptor Accessory Protein antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using IL-18 Receptor Accessory Protein antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human brain using IL-18 Receptor Accessory Protein antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.