

LTB Rabbit pAb

ABZ1556


Product Overview

Name	LTB Rabbit pAb
Accession(Primary)	Q06643
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:200; IF-1:50-1:200; ELISA-1:10000-1:20000


Target information (Q06643)

Gene Name LTB


About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded Human liver using LTB antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using LTB antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human spleen using LTB antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.