

CD66a/b/c Rabbit pAb

ABZ1557


Product Overview

Name	CD66a/b/c Rabbit pAb
Accession(Primary)	P13688
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:50-1:200; ELISA-1:10000-1:20000; IF-1:50-1:200


Target information (P13688)

Gene Name CEACAM1/CEACAM6/CEACAM8


About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded Human kidney using CD66a/b/c antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using CD66a/b/c antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human stomach using CD66a/b/c antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.