


IL-15 Rabbit pAb

ABZ1558


Product Overview

Name	IL-15 Rabbit pAb
Accession(Primary)	P40933
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:50-1:200; ELISA-1:10000-1:20000; IF-1:50-1:200


Target information (P40933)

Gene Name IL15**About Gene Symbol**

Immunohistochemistry analysis of paraffin-embedded Human breast cancer using IL-15 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using IL-15 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human brain using IL-15 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.