

MYH6/7 Rabbit pAb

ABZ1559


Product Overview

Name	MYH6/7 Rabbit pAb
Accession(Primary)	P13533
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:50-1:200; ELISA-1:10000-1:20000; IF-1:50-1:200


Target information (P13533)

Gene Name MYH6/7


About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded Human heart using MYH6/7 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using MYH6/7 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human skeletal muscle using MYH6/7 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.