


B4GALT5 Rabbit pAb


ABZ1560

Product Overview

Name	B4GALT5 Rabbit pAb
Accession(Primary)	O43286
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 45 kDa; Observed MW: 40 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000

Target information (O43286)

Gene Name B4GALT5**About Gene Symbol**



Western blot analysis of B4GALT5 in 293 lysates
using β 1, 4GalT5 antibody.