

LIF Rabbit pAb


ABZ1561

Product Overview


Name	LIF Rabbit pAb
Accession(Primary)	P15018
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 22 kDa; Observed MW: 25 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (P15018)

Gene Name	LIF
About Gene Symbol	



Western blot analysis of LIF in K562 lysates using LIF antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon using LIF antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.