

KLF12 Rabbit pAb


ABZ1562

Product Overview


Name	KLF12 Rabbit pAb
Accession(Primary)	Q9Y4X4
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 44 kDa; Observed MW: 44 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (Q9Y4X4)

Gene Name	KLF12
About Gene Symbol	



Western blot analysis of KLF12 in rat muscle lysates using KLF12 antibody.



Immunohistochemistry analysis of paraffin-embedded mouse brain using KLF12 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.