

FGF5 Rabbit pAb

ABZ1563

Product Overview


Name	FGF5 Rabbit pAb
Accession(Primary)	P12034
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 30 kDa; Observed MW: 30 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (P12034)

Gene Name	FGF5
About Gene Symbol	



Western blot analysis of FGF5 in mouse liver lysates using FGF5 antibody.



Immunohistochemistry analysis of paraffin-embedded Human brain using FGF5 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.