

NRG2 Rabbit pAb


ABZ1564

Product Overview


Name	NRG2 Rabbit pAb
Accession(Primary)	O14511
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 92 kDa; Observed MW: 92 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (O14511)

Gene Name	NRG2
About Gene Symbol	



Western blot analysis of NRG2 in HT-29 lysates using NRG2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human brain using NRG2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.