

Neuropilin 2 Rabbit pAb


ABZ1565

Product Overview


Name	Neuropilin 2 Rabbit pAb
Accession(Primary)	O60462
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 105 kDa; Observed MW: 100 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (O60462)

Gene Name	NRP2
About Gene Symbol	



Western blot analysis of Neuropilin 2 in HepG2 lysates using Neuropilin 2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human brain using Neuropilin 2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.