

CD73 Rabbit pAb


ABZ1566

Product Overview


Name	CD73 Rabbit pAb
Accession(Primary)	P21589
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 63 kDa; Observed MW: 68 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (P21589)

Gene Name	NT5E
About Gene Symbol	



Western blot analysis of CD73 in MCF-7, L929, K562 lysates using CD73 antibody.



Immunohistochemistry analysis of paraffin-embedded Human skin using CD73 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.