

NKG2A/C Rabbit pAb

ABZ1568

Product Overview


Name	NKG2A/C Rabbit pAb
Accession(Primary)	P26715
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 26 kDa; Observed MW: 26 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (P26715)

Gene Name	KLRC1
About Gene Symbol	



Western blot analysis of CD159a/c in L929 lysates using CD159a/c antibody.



Immunohistochemistry analysis of paraffin-embedded Human kidney using CD159a/c antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.