

Placental Alkaline Phosphatase Rabbit pAb


ABZ1569

Product Overview


Name	Placental Alkaline Phosphatase Rabbit pAb
Accession(Primary)	P05187
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 57 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:20000

Target information (P05187)

Gene Name	ALPP/ALPPL2
About Gene Symbol	ALPP (Alkaline Phosphatase, Placental) is a protein phosphatase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Placental Alkaline Phosphatase in 293 lysates using Placental Alkaline Phosphatase antibody.



Western blot analysis of Placental Alkaline Phosphatase in lysates using ALPP/ALPPL2 antibody.

1 2 3 4 customer's cell
6 293T
7 HeLa
8 K562