

FGF16 Rabbit pAb

ABZ1570

Product Overview


Name	FGF16 Rabbit pAb
Accession(Primary)	O43320
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 24 kDa; Observed MW: 24 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (O43320)

Gene Name	FGF16
About Gene Symbol	



Western blot analysis of FGF16 in K562 lysates using FGF16 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung using FGF16 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.