

Histone H2B Rabbit pAb

ABZ1571


Product Overview

Name	Histone H2B Rabbit pAb
Accession(Primary)	P57053
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Cat
Calculated MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:5000


Target information (P57053)

Gene Name	HIST1H2BC/HIST1H2BE/HIST1H2BF/HIST1H2BG/HIST1H2BI/ HIST1H2BD/HIST1H2BH/HIST1H2BK/HIST1H2BL/HIST1H2BM /HIST1H2BN/HIST2H2BF/H2BFS
------------------	---

About Gene Symbol



Western blot analysis of Histone H2B in HeLa lysates using Histone H2B antibody.



Immunohistochemistry analysis of paraffin-embedded Human Breastcancer using Histone H2B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.