

Tenascin C Rabbit pAb


ABZ1573

Product Overview


Name	Tenascin C Rabbit pAb
Accession(Primary)	P24821
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 241 kDa; Observed MW: 241 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:10000

Target information (P24821)

Gene Name	TNC
About Gene Symbol	



Western blot analysis of Tenascin in KB lysates using Tenascin antibody.



Immunohistochemistry analysis of paraffin-embedded Human kidney using Tenascin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.