

CD42d Rabbit pAb


ABZ1574

Product Overview


Name	CD42d Rabbit pAb
Accession(Primary)	P40197
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 61 kDa; Observed MW: 61 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:10000

Target information (P40197)

Gene Name	GP5
About Gene Symbol	



Western blot analysis of CD42d in HeLa lysates using CD42d antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsils using CD42d antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.