

CEACAM6 Rabbit pAb

ABZ1575


Product Overview

Name	CEACAM6 Rabbit pAb
Accession(Primary)	P40198
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 27,37 kDa; Observed MW: 30 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000


Target information (P40198)

Gene Name CEACAM3/CEACAM6

About Gene Symbol



Western blot analysis of CEACAM6 in HepG2, mouse brain lysates using CEACAM6 antibody.



Western blot analysis of CEACAM6 in HEPG2
lysates using CD66c/d antibody.