

CD96 Rabbit pAb

ABZ1576


Product Overview

Name	CD96 Rabbit pAb
Accession(Primary)	P40200
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 66 kDa; Observed MW: 66 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000


Target information (P40200)

Gene Name CD96

About Gene Symbol



Western blot analysis of CD96 in 293UV lysates using CD96 antibody.



Western blot analysis of CD96 in 293UV lysates
using CD96 antibody.