

SIRP alpha/beta Rabbit pAb


ABZ1578

Product Overview


Name	SIRP alpha/beta Rabbit pAb
Accession(Primary)	P78324
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 55 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:10000

Target information (P78324)

Gene Name	SIRPA/SIRPB1
About Gene Symbol	



Western blot analysis of SIRP alpha/beta in HeLa lysates using SIRP alpha/beta antibody.



Immunohistochemistry analysis of paraffin-embedded Human brain using SIRP alpha/beta antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.