

CD232 Rabbit pAb

ABZ1579


Product Overview

Name	CD232 Rabbit pAb
Accession(Primary)	O60486
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 176 kDa; Observed MW: 176 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000


Target information (O60486)

Gene Name PLXNC1**About Gene Symbol**

3T3



Western blot analysis of CD232 in 3T3 lysates using CD232 antibody.



Western blot analysis of CD232 in 3T3 lysates
using CD232 antibody.