

Poliovirus Receptor/PVR Rabbit pAb


ABZ1580

Product Overview


Name	Poliovirus Receptor/PVR Rabbit pAb
Accession(Primary)	P15151
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 45 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:500-1:2000; IF-1:100-1:300

Target information (P15151)

Gene Name	PVR
About Gene Symbol	PVR (PVR Cell Adhesion Molecule) is a membrane receptor involved in adaptive immunity. It is an important research target in infectious disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of CD155 in KB lysates using CD155 antibody.



Western blot analysis of CD155 in KB lysates using CD155 antibody.