

CD179b Rabbit pAb


ABZ1581

Product Overview


Name	CD179b Rabbit pAb
Accession(Primary)	P15814
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 23 kDa; Observed MW: 23 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:10000

Target information (P15814)

Gene Name	IGLL1
About Gene Symbol	



Western blot analysis of CD179b in KB lysates using CD179b antibody.



Immunohistochemistry analysis of paraffin-embedded Human lymph using CD179b antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.