

LEFTY Rabbit pAb

ABZ1582


Product Overview


Name	LEFTY Rabbit pAb
Accession(Primary)	O75610
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300

Target information (O75610)

Gene Name LEFTY1/LEFTY2

About Gene Symbol





Western blot analysis of LEFTY in HEPG2UV
lysates using Lefty antibody.