


Epigen Rabbit pAb

ABZ1583

Product Overview

Name	Epigen Rabbit pAb
Accession(Primary)	Q6UW88
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 17 kDa; Observed MW: 17 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000

Target information (Q6UW88)

Gene Name EPGN**About Gene Symbol**

Western blot analysis of Epigen in 293T, 823 lysates using Epigen antibody.

