

MOB3A/B Rabbit pAb


ABZ1584

Product Overview


Name	MOB3A/B Rabbit pAb
Accession(Primary)	Q96BX8
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 25 kDa; Observed MW: 25 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:10000

Target information (Q96BX8)

Gene Name	MOB3A
About Gene Symbol	



Western blot analysis of MOB3A/B in 3T3 lysates using MOB3A/B antibody.



Immunohistochemistry analysis of paraffin-embedded Human skin using MOB3A/B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.