


MRPL35 Rabbit pAb

ABZ1585


Product Overview

Name	MRPL35 Rabbit pAb
Accession(Primary)	Q9NZE8
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 22 kDa; Observed MW: 22 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300

Target information (Q9NZE8)

Gene Name MRPL35**About Gene Symbol**

Western blot analysis of MRPL35 in 293T , COLO
lysates using MRPL35 antibody.



Western blot analysis of MRPL35 in 293 lysates
using MRPL35 antibody.