

Factor XIII Rabbit pAb


ABZ1586

Product Overview


Name	Factor XIII Rabbit pAb
Accession(Primary)	P05160
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 76 kDa; Observed MW: 65 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:10000

Target information (P05160)

Gene Name	F13B
About Gene Symbol	



Western blot analysis of Factor XIII in various lysates using Factor XIII antibody.



Immunohistochemistry analysis of paraffin-embedded Human liver using Factor XIII antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.