


Integrin alpha 3 Rabbit pAb

ABZ1587


Product Overview

Name	Integrin alpha 3 Rabbit pAb
Accession(Primary)	P26006
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 117 kDa; Observed MW: 117 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:40000; IHC-1:50-1:100

Target information (P26006)

Gene Name ITGA3**About Gene Symbol**

Western blot analysis of Integrin alpha 3 in hela lysates using Integrin alpha 3 antibody.



Immunohistochemistry analysis of paraffin-embedded mouse spleen using Integrin alpha 3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.