

SOCS3 Rabbit pAb

ABZ1588


Product Overview

Name	SOCS3 Rabbit pAb
Accession(Primary)	O14543
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	WB-1:500-1:2000; IHC:1:50-1:300; IF-1:50-1:200


Target information (O14543)

Gene Name SOCS3

About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using SOCS3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using SOCS3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.