

## CCR1 Rabbit pAb

ABZ1589


### Product Overview

|                             |                                                         |
|-----------------------------|---------------------------------------------------------|
| <b>Name</b>                 | CCR1 Rabbit pAb                                         |
| <b>Accession(Primary)</b>   | P32246                                                  |
| <b>Host Species</b>         | Rabbit                                                  |
| <b>Clonality</b>            | Polyclonal                                              |
| <b>Isotype</b>              | IgG                                                     |
| <b>Species Reactivity</b>   | Human, Mouse                                            |
| <b>Recommended Dilution</b> | IHC-1:50-1:200; ELISA-1:10000-1:20000;<br>IF-1:50-1:200 |


### Target information (P32246)

**Gene Name** CCR1

**About Gene Symbol**



Immunohistochemistry analysis of paraffin-embedded Human small intestine using CCR1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using CKR 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.