


CD300b Rabbit pAb

ABZ1590


Product Overview

Name	CD300b Rabbit pAb
Accession(Primary)	A8K4G0
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Recommended Dilution	IHC-1:50-1:200; ELISA-1:10000-1:20000; IF-1:50-1:200

Target information (A8K4G0)

Gene Name CD300LB**About Gene Symbol**

Immunohistochemistry analysis of paraffin-embedded Human liver cancer using CD300b antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using CD300b antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.