


NNT1 Rabbit pAb

ABZ1591


Product Overview

Name	NNT1 Rabbit pAb
Accession(Primary)	Q9UBD9
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:50-1:200; ELISA-1:10000-1:20000; IF-1:50-1:200

Target information (Q9UBD9)

Gene Name CLCF1**About Gene Symbol**

Immunohistochemistry analysis of paraffin-embedded Human brain using NNT1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using NNT1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.