

CDNF Rabbit pAb

ABZ1592


Product Overview

Name	CDNF Rabbit pAb
Accession(Primary)	Q49AH0
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:50-1:200; ELISA-1:10000-1:20000; IF-1:50-1:200


Target information (Q49AH0)

Gene Name CDNF

About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded Human liver using CDNF antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using CDNF antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.