

5HT3A receptor/5ht3 Rabbit pAb


ABZ1597

Product Overview

Name	5HT3A receptor/5ht3 Rabbit pAb
Accession(Primary)	P46098
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 55 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:10000

Target information (P46098)

Gene Name	HTR3A
About Gene Symbol	



Western blot analysis of HTR3A in various lysates using HTR3A antibody.