

Ubiquitin Rabbit pAb


ABZ1598

Product Overview

Name	Ubiquitin Rabbit pAb
Accession(Primary)	P0CG47
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Refer to figures
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (P0CG47)

Gene Name	UBB
About Gene Symbol	



Western blot analysis of Ubiquitin in various lysates using Ubiquitin antibody.