

ALPK1 Rabbit pAb


ABZ1599

Product Overview

Name	ALPK1 Rabbit pAb
Accession(Primary)	Q96QP1
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 139 kDa; Observed MW: 139 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:20000

Target information (Q96QP1)

Gene Name	ALPK1
About Gene Symbol	



Western blot analysis of ALPK1 in various lysates
using ALPK1 antibody.